

DATASHEET

The Forceps/Vacuum Delivery OB Manikin has all of the same realistic and accurate anatomical manikin features as the Obstetrical Manikin (110-180U), but with added features of premature and full term fetuses.

The soft vinyl pelvis replicates the resistance encountered in a delivery requiring forceps or vacuum intervention.

The included abdominal overlay is different from those that come with the Obstetrical Manikin, and can be interchanged or removed.

A forceps or vacuum extractor is not included.

Forceps/Vacuum Delivery OB Manikin



KEY BENEFITS

- Same features as the Obstetrical Manikin with focus on forceps and vacuum delivery
- Realistic pelvic cavity
- Interchangeable abdominal covers
- Soft vinyl pelvis replicates resistance during labor

ORDERING SPECIFICATIONS

Item no.

110-1377U

Forceps/Vacuum Delivery OB Manikin

incl. abdominal overlay, extra vulva, blood powder

Ambu

🏢 Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu.com